

Technical Data Sheet Issue Date: 01-01-03		Product Code: VT-116 Glue Remover
--	---	--

Product Description:

“VitalFix” **Glue Remover** is an ideal solution to solve most common problems associated with careless use of glue.

Applications:

It removes Super Glue from accidental bonded finger. Other applications include removing stained surfaces, fabrics, cloth and badly bonded objects. It can be used to detach also traces of sticky labels, inks and permanent marker stained.

Properties:

Chemical Type	Industrial Debonder
Solvent	Butyrolactone
Boiling Point	204°C
Appearance	Clear/Amber Liquid
Flashpoint	93°C
Density	1, 12

Directions:

1. Puncture with top of cap and install nozzle.
2. **Stained Surfaces/Stained Clothes/Sticky Labels/Badly Bonded Objects**
Apply the solution onto surface and place a folded tissue on top. Remove tissue and wash with water. Repeat operation if require.
3. **Bonded Finger**
Apply solution onto bonded area and allow several minutes to let glue soften. Gently peel bonded area apart allowing the solution to penetrate and break down the bond. Do not attempt to force the skin apart by pulling. Carefully separate the fingers and rinse with water.
4. Replace cap after use and store in cool, dry place.

CAUTION:

Contains gamma-Butyrolactone. Irritating to eyes. Harmful if swallowed. In case of skin contact, soak in warm soapy water, remove contaminated clothing immediately. In case of eye contact, flush with cold water for 15 minutes and seek medical attention immediately.

Storage:

Keep out of reach of children. Use in well ventilated areas.

Every endeavor has been made to ensure that the information given herein is true and reliable but it is given only for the guidance of our customers. The company cannot accept any responsibility for the loss or damage that may result from the use of the information, due to the possibility of various of processing or working conditions and of workmanship outside our control. Users are advised to confirm suitability of this product by their own tests.

- End -